



NCSA Technology Procurement 2023-2024

Effective March 16, 2023 - March 15, 2024

Zoll Medical Corporation

Medical Devices and Technology Solutions for Law Enforcement: Catalog Offering

8% Discount off MSRP for Zoll Medical Equipment

Key Notes:

- Freight & Sales Tax are not included in pricing
- Services are available please contact Paula Sauls or Tony Bailey for details
- MSRP prices are subject to change

The list below represents popular models. Please contact ProLogic ITS if the item you seek is not listed.

AUTOMATED EXTERNAL DEFIBRILLATORS

21000010102011010	ZOLL AED PLUS, PS SERIES, LCD, NO VOICE RCDG,
	GEMS IT, #2012511-001, ENGLISH
21400010101011010	AED PLUS, PS SERIES, W/BASIC CVR, LCD, NO VOICE
	RCDG, 60HZ, ENG
21400010102011010	AED PLUS, PS SERIES, W/AED CVR, LCD, NO VOICE
	RCDG, 60HZ, ENG
21400010302011010	AED PLUS, PS SERIES NO GRAPHICS, W/AED CVR,
	LCD, NO VOICE RCDG, 60HZ, ENG
22200000102011010	AED PLUS, PS SERIES, W/AED CVR, LCD, NO VOICE
	RCDG, POLICE
22200010101011010	AED PLUS, PS SERIES, W/BASIC CVR, LCD, NO VOICE
	RCDG, 60HZ, ENGLISH, PLUSRX
22200010102011010	AED PLUS, PS SERIES, W/AED CVR, LCD, NO VOICE
	RCDG, 60HZ, ENGLISH, PLUSRX
22200010302011010	AED PLUS, PS SERIES NO GRAPHICS, W/AED CVR,
	LCD, NO VOICE RCDG, 60HX, ENGLISH, PLUSRX
2319-001-Z	VIRTUAL CPR TRAINING
MOBILIZE RESCUE SYSTEMS	
8911-001000-01	MOBILIZE RESCUE SYSTEMS, COMPREHENSIVE

For quotes or more information, contact:

Paula Sauls, Account Manager – <u>Paula.Sauls@ProLogicITS.com</u> Tony Bailey, Government Manager – <u>Tony.Bailey@ProLogicITS.com</u>





NCSA Technology Procurement 2023-2024

Effective March 16, 2023 - March 15, 2024

8911-002000-01	MOBILIZE RESCUE SYSTEMS, MOBILE
8911-003000-01	MOBILIZE RESCUE SYSTEMS, COMPACT
8911-003500-01	MOBILIZE RESCUE SYSTEMS, TRAINING KIT
8911-004000-01	MOBILIZE RESCUE SYSTEMS, UTILITY
8911-005000-01	MOBILIZE RESCUE SYSTEMS, PUBLIC ACCESS RESCUE STATION
90110200499991010	AED PRO, SEMI AUTO W/MNL OVERIDE, NO BATTERY, NO ELECTRODES
90110400499991010	AED PRO, SEMI AUTOMATIC ONLY, NO BATTERY, NO ELECTRODES
90110600499991010	AED PRO, MANUAL ONLY, NO VOICE, LCD, ENGLISH
RESQ SYSTEM	
12-0823-000	RESQPUMP ACD-CPD DEVICE
12-0825-000	RESQCPR SYSTEM


For quotes or more information, contact:

Paula Sauls, Account Manager – Paula.Sauls@ProLogicITS.com Tony Bailey, Government Manager – <u>Tony.Bailey@ProLogicITS.com</u>



Survival increases with early intervention

Research shows that the probability of survival increases dramatically when CPR is performed and when an AED is applied before an ambulance arrives.¹



Probability of Survival (Discharged from hospital and alive 30 days later)

Where is the AED?

Too often, the answer is, "We don't have one." Sadly, the same research that demonstrated a nearly fivefold increase in survivability (from 5% to 24%) when an AED is used, also showed that an AED is available only 2% of the time.

Prol\sigic ITS POWERHEART G5



Automatic Self-Testing

battery and internal circuitry

to ensure your AED is always

Rescue Ready®.

and recoil.

Conducts daily self-tests on pads,

Intellisense™ CPR Feedback

CPR prompts for rate, depth

Four distinct text and voice corrective

Intellisense[™] Defibrillation Pad

Non-polarized pads can be placed in either location on the chest for simple, fast placement.

RescueCoach™

User-paced voice and text prompts guide rescuers through each critical step of the rescue, including CPR instruction.

USB Data Download

90 minutes of rescue data available via USB download.

Dual Language

Change from primary to secondary language with the push of a button, without delaying or interrupting the rescue.


Whether you are an experienced rescuer or a first-time responder, you will appreciate how the Powerheart G5 AED combines ease of use, reliability, and advanced technologies to help victims survive sudden cardiac arrest.

The Powerheart G5 AED with Intellisense CPR provides real-time, Guideline-driven feedback to ensure delivery of high-quality compressions. One-button dual-language functionality and daily self-tests give rescuers the confidence to act.

ZOLL®

POWERHEART 65







DESIGNED FOR UNEXPECTED HEROES™

The survival rate from sudden cardiac arrest (SCA) is less than 1% worldwide.⁴ According to the American Heart Association (AHA), high-quality CPR, along with early defibrillation, can help improve SCA survival rates.⁵

The Powerheart® G5 AED is the only AED to combine real-time CPR feedback with dual-language functionality, user-paced rescue prompting, and fast shock times.



-66

"I was not expecting the AED to walk me through doing compressions. I was sure thankful for that 'press, press, press'."

Mark Tucker, Rescuer CM Tucker Lumber in Pageland, SC

•

It happens in a split second: A person collapses – the victim of sudden cardiac arrest – and the clock starts ticking. You call for help and medical personnel are on their way. But precious seconds are racing by. The victim's best chance for survival depends on immediate CPR and a rescuer with an automated external defibrillator (AED).



CPR plus early defibrillation can **more than double the** rate of survival from out-of-hospital cardiac arrest.

-2015 AHA Guidelines²

THE POWER TO HELP

The latest guidelines from resuscitation councils worldwide are clear: Successful defibrillation must be supported by high-quality CPR.

The AED's first heart analysis is critical. If it calls for no shock, only high-quality CPR can lead to a shockable rhythm on the next heart analysis. High-quality CPR provides a struggling heart with oxygenated blood needed to help restore a normal rhythm.



"...monitoring of CPR quality is arguably one of the most significant advances in resuscitation practice in the past 20 years and should be incorporated into every resuscitation and every professional rescuer program."

- American Heart Association CPR Quality Consensus Statement³



Real-Time CPR Feedback

The Powerheart G5 AED uses Intellisense CPR (ICPR) technology to monitor CPR chest compressions, providing rescuers with four distinct corrective voice and text prompts: Press Slower, Press Faster, Press Softer, and Press Harder and Fully Release.



User-Paced Rescue Prompts

RescueCoach® prompts guide users step-by-step through the rescue process. The AED recognizes actions taken and ensures that the rescuer completes every critical task before moving on to the next.




Rescue Ready®

The Powerheart G5 AED performs comprehensive daily self-testing to ensure it is ready to rescue. Studies have shown that approximately 46% of AED failures are due to pad, pad connector, or battery power problems. The Powerheart AED self-tests detect and help prevent these issues.




Dual-Language Capabilities

With a touch of a button, CPR Instruction, ICPR corrections, and AED prompts can quickly switch from the primary language to a pre-programmed alternate language without delaying the rescue.



99



Your officers are often the first to arrive to a medical emergency. Some emergencies, like severe bleeding and airway complications, can become fatal within a few minutes — before EMS arrives. Do your officers have the resources and knowledge to help at the scene? ZOLL® can help make sure they do.

Mobilize[™] from ZOLL is a set of interactive trauma and medical emergency kits designed to empower bystanders to help victims during unexpected medical emergencies.



Mobilize Rescue app: Helps untrained bystanders provide high-quality medical care by following the real-time diagnostic prompts in the interactive app. The algorithm and treatment directions are based on protocols and best practices from the American Heart Association, Tactical Emergency Casualty Care (TECC), and the American College of Surgeons.



High-quality supplies: Mobilize trauma kits contain modern medical supplies to treat and manage a wide range of emergencies from severe bleeding and chest trauma to seizures, allergic reactions, and more.

¹National Trauma Institute website, www.nattrauma.org/what-is-trauma/trauma-statistics-fact/. Accessed June 24, 2020 ²U.S. Department of Homeland Security, www.dhs.gov/stopthebleed ³https://www.nattrauma.org/what-is-trauma/trauma-statistics-facts/



Trauma is the number one cause of death in ages 1 to 46, accounting for

30% OF ALL LIVES LOST IN THE US¹



Blood loss must be stopped ASAP. Bleeding can become fatal in as little as

3-5 MINUTES²



30,000
TRAUMATIC DEATHS

could be prevented each year with increased access to trauma care.³

Phone: (866) 923-0513 Web: www.prologicits.com

Product Listings:

Mobilize Comprehensive trauma system

Hardened, water-resistant carry case with integrated Mobilize Rescue app and complete set of rescue supplies. (8911-001000-01)

Part #

Sale Price



Mobilize Mobile trauma kit

Includes the same contents as the Comprehensive rescue system in a durable military-grade backpack. (8911-002000-01)

Part #

Sale Price



Mobilize Public Access trauma kit

Four Utility kits* in a soft carry case with unlimited app licenses. Each Utility kit includes supplies for trauma, bleeding control, CPR, and more. (8911-005000-01)

Part #

Sale Price

*Individual replacement kits available



ProL<u>§gic ITS</u>

Mobilize **Compact** bleeding control and airway management kit

Includes multi-user app license and supplies for trauma, bleeding control, CPR, and more in a durable, military-grade case. (8911-003000-01)

Part #

Sale Price

Cabinets:

Wall Cabinet for Public Access

Alarmed: 8911-000510-01

Non-alarmed: 8911-000500-01

Alarmed Wall Cabinet for Comprehensive or Mobile

8911-000470-01



Want to know more?

Please contact your ProLogic ITS Law Enforcement Market Specialist:

